

# Disposable SpO<sub>2</sub> Sensors



## Features

- *Engineered to function equivalent to OEMs*
- *Quality approved with CE Certificate and FDA 510K*
- *Single patient use eliminates cross-contamination issues*
- *Reliable performance with sensitive and accurate responses*
- *Cost effective*

Available across the entire line of primary brand compatible disposable sensors

- |                    |                       |                |                 |
|--------------------|-----------------------|----------------|-----------------|
| ▪ <b>BCI/Smith</b> | ▪ <b>CSI</b>          | ▪ <b>GE</b>    | ▪ <b>Masimo</b> |
| ▪ <b>Nellcor</b>   | ▪ <b>Nihon Kohden</b> | ▪ <b>Nonin</b> |                 |



PN: SA-001X-XC  
Price: \$20.00/pc


## Oximax Adapter Cable

The Oximax Adapter Cable cuts the cost of Oximax disposable sensors allowing you to use the non-Oximax Nellcor disposable sensors with same performance functionally as the Oximax sensors.


## Nellcor- Oximax

Cable length: 0.9m

Package: 24 pcs/box



SSD-001XC-W09A



SSD-001XC-S09PI



Sensor Connector

| Part Number   | SSD-001XC-W09A           | SSD-001XC-W09P       | SSD-001XC-W09I   | SSD-001XC-W09AN                | SSD-001XC-S09AN                | SSD-001XC-S09PI            |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM#          | MAXA, MAXAL              | MAXP                 | MAXI             | MAXAN                          | MAXA, MAXI                     | MAXP, MAXI                 |
| Patient Sizes | Adult<br>>40kg           | Pediatric<br>10-50kg | Infant<br>3-20kg | Neonate/Adult<br><3kg or >40kg | Adult/Neonate<br>>40kg or <3kg | Pediatric/Infant<br>3-40kg |
| Tapes         | Adhesive Non-woven fiber |                      |                  |                                | Non-adhesive foam              |                            |

## Nellcor- Oxismart (non Oximax)

Cable length: 0.9m


Package: 24 pcs/box



SSD-001-W09A



SSD-001-S09AN



Sensor Connector

| Part Number   | SSD-001-W09A             | SSD-001-W09P         | SSD-001-W09I     | SSD-001-W09AN                  | SSD-001-S09AN                  | SSD-001-S09PI              |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM#          | D25                      | D20                  | I20              | N25                            | N25, D25                       | D20, I20                   |
| Patient Sizes | Adult<br>>40kg           | Pediatric<br>10-50kg | Infant<br>3-20kg | Neonate/Adult<br><3kg or >40kg | Adult/Neonate<br>>40kg or <3kg | Pediatric/Infant<br>3-40kg |
| Tapes         | Adhesive Non-woven fiber |                      |                  |                                | Non-adhesive foam              |                            |

## Masimo- LNCS

Cable length: 0.9m

Package: 20 pcs/box



SSD-005RB-W09A

SSD-005RB-S09AN

Sensor Connector

| Part Number   | SSD-005RB-W09A           | SSD-005RB-W09P       | SSD-005RB-W09I   | SSD-005RB-W09AN                | SSD-005RB-S09AN                | SSD-005RB-S09PI            |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM#          | 1859,2317                | 1860,2318            | 1861,2328        | 2329,1862                      | 2606                           | 1860,2328                  |
| Patient Sizes | Adult<br>>40kg           | Pediatric<br>10-50kg | Infant<br>3-20kg | Neonate/Adult<br><3kg or >40kg | Adult/Neonate<br>>40kg or <3kg | Pediatric/Infant<br>3-40kg |
| Tapes         | Adhesive Non-woven fiber |                      |                  |                                | Non-adhesive foam              |                            |

## Masimo- M-LNCS

Cable length: 0.9m

Package: 20 pcs/box



SSD-007RB-W09A

SSD-007RB-S09AN

Sensor Connector

| Part Number   | SSD-007RB-W09A           | SSD-007RB-W09P       | SSD-007RB-W09I   | SSD-007RB-W09AN                | SSD-007RB-S09AN                | SSD-007RB-S09PI            |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM#          | 2508,2509                | 2510,2511            | 2512,2513        | 2514,2515                      | 2508,2509                      | 2510,2511                  |
| Patient Sizes | Adult<br>>40kg           | Pediatric<br>10-50kg | Infant<br>3-20kg | Neonate/Adult<br><3kg or >40kg | Adult/Neonate<br>>40kg or <3kg | Pediatric/Infant<br>3-40kg |
| Tapes         | Adhesive Non-woven fiber |                      |                  |                                | Non-adhesive foam              |                            |

## Masimo- LNOP

Cable length: 0.9m

Package: 20 pcs/box



SSD-006RB-W09A

SSD-006RB-S09AN


Sensor Connector

| Part Number   | SSD-006RB-W09A           | SSD-006RB-W09P       | SSD-006RB-W09I   | SSD-006RB-W09AN                | SSD-006RB-S09AN                | SSD-006RB-S09PI            |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM#          | 1001                     | 1025                 | 1800             | 2358                           | 1001                           | 1025                       |
| Patient Sizes | Adult<br>>40kg           | Pediatric<br>10-50kg | Infant<br>3-20kg | Neonate/Adult<br><3kg or >40kg | Adult/Neonate<br>>40kg or <3kg | Pediatric/Infant<br>3-40kg |
| Tapes         | Adhesive Non-woven fiber |                      |                  |                                | Non-adhesive foam              |                            |

## GE- Datex

Cable length: 0.9m

Package: 24 pcs/box



SSD-013-W09A

SSD-013-S09AN


Sensor Connector

| Part Number   | SSD-013-W09A             | SSD-013 -W09P        | SSD-013-W09I     | SSD-013-W09AN                  | SSD-013-S09AN                  | SSD-013-S09PI              |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM#          | SAS-AF                   | SAS-AP               | SAS-AP           | SAS-AF                         | SAS-AF                         | SAS-AP                     |
| Patient Sizes | Adult<br>>40kg           | Pediatric<br>10-50kg | Infant<br>3-20kg | Neonate/Adult<br><3kg or >40kg | Adult/Neonate<br>>40kg or <3kg | Pediatric/Infant<br>3-40kg |
| Tapes         | Adhesive Non-woven fiber |                      |                  |                                | Non-adhesive foam              |                            |

## GE- TruSignal

Cable length: 0.9m

Package: 24 pcs/box



SSD-018-W09A

SSD-018-S09AN

Sensor Connector

| Part Number   | SSD-018-W09A             | SSD-018-W09P         | SSD-018-W09I     | SSD-018-W09AN                  | SSD-018-S09AN                  | SSD-018-S09PI              |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM#          | TS-AF, TS-AAW            | TS-AP, TS-PAW        | TS-AF            | TS-AF                          | TS-AF, TS-AAW                  | TS-AP, TS-PAW              |
| Patient Sizes | Adult<br>>40kg           | Pediatric<br>10-50kg | Infant<br>3-20kg | Neonate/Adult<br><3kg or >40kg | Adult/Neonate<br>>40kg or <3kg | Pediatric/Infant<br>3-40kg |
| Tapes         | Adhesive Non-woven fiber |                      |                  |                                | Non-adhesive foam              |                            |

## GE- OxyTip

Cable length: 0.9m

Package: 24 pcs/box



SSD-019-W09A

SSD-019-S09AN

Sensor Connector

| Part Number   | SSD-019-W09A             | SSD-019-W09P         | SSD-019-W09I     | SSD-019-W09AN                  | SSD-019-S09AN                  | SSD-019-S09PI              |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM#          | OXY-AF                   | OXY-AP,OXY-AF        | OXY-AF           | OXY-AF                         | OXY-AF                         | OXY-AP,OXY-AF              |
| Patient Sizes | Adult<br>>40kg           | Pediatric<br>10-50kg | Infant<br>3-20kg | Neonate/Adult<br><3kg or >40kg | Adult/Neonate<br>>40kg or <3kg | Pediatric/Infant<br>3-40kg |
| Tapes         | Adhesive Non-woven fiber |                      |                  |                                | Non-adhesive foam              |                            |

## ■ Nihon Kohden

Cable length: 0.9m


Package: 24 pcs/box



SSD-043-W09A



SSD-043-S09AN



Sensor Connector

| Part Number   | SSD-043-W09A             | SSD-043-W09P         | SSD-043-W09I     | SSD-043-W09AN                  | SSD-043-S09AN                  | SSD-043-S09PI              |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM#          | TL-271                   | TL-272               | TL-274           | TL-273                         | TL-273                         | TL-272/T274                |
| Patient Sizes | Adult<br>>40kg           | Pediatric<br>10-50kg | Infant<br>3-20kg | Neonate/Adult<br><3kg or >40kg | Adult/Neonate<br>>40kg or <3kg | Pediatric/Infant<br>3-40kg |
| Tapes         | Adhesive Non-woven fiber |                      |                  |                                | Non-adhesive foam              |                            |

## ■ Nonin

Cable length: 0.9m


Package: 24 pcs/box



SSD-040-W09A



SSD-040-S09PI




Sensor Connector

| Part Number   | SSD-040-W09A             | SSD-040-W09P         | SSD-040-W09I     | SSD-040-W09AN                  | SSD-040-S09AN                  | SSD-040-S09PI              |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM#          | 7000A, 6000CA            | 7000P, 6000CP        | 7000I, 6000CI    | 7000N, 6000CN                  | 7000A, 7000N,                  | 7000P, 7000I,              |
| Patient Sizes | Adult<br>>40kg           | Pediatric<br>10-50kg | Infant<br>3-20kg | Neonate/Adult<br><3kg or >40kg | Adult/Neonate<br>>40kg or <3kg | Pediatric/Infant<br>3-40kg |
| Tapes         | Adhesive Non-woven fiber |                      |                  |                                | Non-adhesive foam              |                            |

## BCI/Smith

Cable length: 0.9m

Package: 24 pcs/box



SSD-036-W09A

SSD-036-S09AN


Sensor Connector

| Part Number   | SSD-036-W09A             | SSD-036-W09P         | SSD-036-W09I     | SSD-036-W09AN                  | SSD-036-S09AN                  | SSD-036-S09PI              |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM#          | 1300                     | 1301                 | 1303             | 1302                           | 1300, 1302                     | 1301, 1303                 |
| Patient Sizes | Adult<br>>40kg           | Pediatric<br>10-50kg | Infant<br>3-20kg | Neonate/Adult<br><3kg or >40kg | Adult/Neonate<br>>40kg or <3kg | Pediatric/Infant<br>3-40kg |
| Tapes         | Adhesive Non-woven fiber |                      |                  |                                | Non-adhesive foam              |                            |

## CSI (Criticare)

Cable length: 0.9m

Package: 24 pcs/box



Sensor Connector

| Part Number   | SSD-039-W09A             | SSD-039-W09P         | SSD-039-W09I     | SSD-039-W09AN                  | SSD-039-S09AN                  | SSD-039-S09PI              |
|---------------|--------------------------|----------------------|------------------|--------------------------------|--------------------------------|----------------------------|
| OEM#          | 570SD                    | 571SD                | 572SD            | 573SD                          | 570SD                          | 571SD, 572SD               |
| Patient Sizes | Adult<br>>40kg           | Pediatric<br>10-50kg | Infant<br>3-20kg | Neonate/Adult<br><3kg or >40kg | Adult/Neonate<br>>40kg or <3kg | Pediatric/Infant<br>3-40kg |
| Tapes         | Adhesive Non-woven fiber |                      |                  |                                | Non-adhesive foam              |                            |