

EKG Patient Cables

- Engineered to functionally equivalent to OEMs
- Quality approved with CE Certificated and FDA 510K
- Being connected the EKG machines and adaptors
- Affordable Prices
- 12 months Warranty



GE Marquette



EKT-GE-CAM14 / EKT-GE-CAM14-S



EKT-GE-AM4 / EKT-GE-AM4-S



EKT-GE-MAXP / EKT-GE-MAXP-S



EKT-GE-MAX1

Part Number	Equipment Models	Acquisition Module	Cable Length	OEM Part Number	Price(USD)
EKT-GE-CAM14	CASE 8000 Exercise testing System	CAM 14	4.6m+0.3m	2016560-002	
EKT-GE-CAM14-S	MAC 5000 Resting System		1.3m+0.3m	2016560-001	
EKT-GE-AM4	CASE 12, CASE 15, MAC PC, MAC 6, MAC 12, MAC 15, CENTRA	AM4 or AM5	4.6m+0.3m	700044-202	
EKT-GE-AM4-S			1.3m+0.3m	700044-203	
EKT-GE-MAXP	CASE 16, MAX-Personal, MAC VU System	AM4 or AM5	4.6m+0.3m	700044-204	
EKT-GE-MAXP-S			1.3m+0.3m	700044-205	
EKT-GE-MAX1	MAX-1, MAX-4, MAX-6	AM4 or AM5	4.6m+0.3m	700044-201	

Orantech Inc.

T +86-755 2369 9939 F +86-755 2324 3594 E info@orantech.com W www.orantech.com

