

Medical Oxygen Sensors



Features

- Engineered to function equivalent to OEMs
- CE marked products
- Precise measurements
- Long service life
- Made in German



Compatible with these major brands

- | | |
|--------------------------------|-----------------------|
| ▪ Analytical Industries | ▪ Datex-Ohmeda |
| ▪ City Technology | ▪ Envitec |
| ▪ Drager | ▪ Hudson-RCI |
| ▪ Hamilton | ▪ Teledyne |
| ▪ Maxtec | |

Medical Oxygen Sensor

Descriptions	OEM Machines & Sensors	Price/USD
 OS-01 Response Time:17s Warranty: 12 months	Analytical Industries(PSR-11-915) City Technologies(MOX-6) Draeger Cells(6803290) Envitec Cells(OOM101) GE-Datex-Ohmeda Cells(0237-2034-700) GE-DatexOhmeda Instrum(Excel 210-7800, Excel 7800, Modulus 7800, 4700 Oxicap, 5100, 5120, 5125, 5140, 5200, 5250, 5250 RMG, 5400, 7200, 7800, 7810 Vent, S5172802) Hamilton Cells(HM-10) Hudson-RCI Cells(5556) Maxtec(MAX-1) MSA Cells(655264, 655263) Nuova(E-15/0) Teledyne Cells(R-15) Teledyne Instruments(#C-41131) Viamed(R-15)	
 OS-02 Response Time:17s Warranty: 12 months	Analytical Industries(PSR-11-915-2) Draeger Cells(6850645, 6803290) Envitec Cells(OOM101) Draeger Instruments(Narkomed 1A,2,2B,3, 4, 6000 Babylog 8000, SA2 Anaesthesia Mach, Evita, Evita 2, Evita 4,Evita 2 Dura Fabius, Fabius GS, Oxydig, Oxycom 100 D PM8030, I-8000,Incubator 8000IC, Anemone) Envitec Cells(OOM201) GE-Datex-Ohmeda Cells(6051-000-224) Hamilton Cells(HM-14, HM-13) Hudson-RCI Cells(5557, V-14) Maxtec(MAX-11) MSA Cells(655263,655264) Nuova(E-15/2) Teledyne Cells(R-23MED) Teledyne Instruments(#C-61339) Viamed(R-23V)	

Descriptions	OEM Machines & Sensors	Price/USD
 OS-03 Response Time: <12s Warranty: 16 months	Envitec Cells(OOM202-1) GE Datex-Engström(Elvira ,Erica) GE-Datex-Ohmeda Cells(6051-0000-217) Hamilton Cells(HM-11) Hudson-RCI Cells(5567, V-11) Maxtec(MAX-13) MSA Cells(406931) MSA Instruments(Mini-Ox I - III) Nuova(E-17/J) Taema Instruments(Eole, Cesar) Viamed(R-41V)	
 OS-04 Response Time: <12s Warranty: 16 months	City Technologies(MOX-9) Envitec Cells(OOM202) GE-Datex-Ohmeda Cells(6051-0000-222) Maxtec(MAX-12) MSA Cells(472062 , 470644) Nuova(E-17/M) Siemens Cells(9004979) Siemens Instruments(Ventilators:Servo 900 /C/D, Servo 700/710, Servo 300, V710 300) Taema Instruments(YR049700, Cesar2, Horus) Teledyne Cells(R-24Med) Teledyne Instruments(#B 62321) Viamed(R-24Med)	
 OS-05 Response Time: <20s Warranty: 12months	Analytical Industries(PSR-11-33) Envitec Cells(OOM107) Hamilton Cells(HM-01) Hamilton Instruments(Veolar, Amadeus, Aladin) Hudson-RCI Cells(5568, C1, V-01) Hudson-RCI Instruments(5525, 5528, 5529, 5550, 5552) Maxtec(MAX-19) MSA Cells(478841,485955) Nuova(E-33/0) Taema Instruments(KV050501, Oxy 2100) Teledyne Cells(T1/T2) Teledyne Instruments(TED 160, TED 190) Viamed(T-1/T-2)	

Descriptions	OEM Machines & Sensors	Price/USD
 OS-06 Response Time: <20s Warranty: 12months	Analytical Industries(PSR-11-33-2) Envitec Cells(OOM107-2) Hudson-RCI Cells(5569, V-02) Maxtec(MAX-18) MSA Cells(485905, 655262) Nuova(E-33/1) Teledyne Cells(C2R) Viamed(C-2R)	
 OS-07 Response Time: <8s Warranty: 16months	Analytical Industries(PSR-11-77) City Technologies(MOX-4) Envitec Cells(OOM105) GE-Datex-Ohmeda Cells(6051-0000-219) Hamilton Cells(HM-03) Hudson-RCI Cells(5558, V-03) Maxtec(MAX-17) MSA Cells(804674) Nuova(E-77/0) Teledyne Cells(T-7) Teledyne Instruments(TED 60T, TED 191T, TED 200T, TED 200T7) Viamed(T-7)	
 OS-07S Response Time: <15s Warranty: 12months	Draeger Instruments(I-8000) Envitec Cells(OOM104) Maxtec(MAX-11i) Nuova(E-8000)	
 OS-08 Response Time: <12s Warranty: 16months	Draeger Instruments(I-8000) Envitec Cells(OOM104) Maxtec(MAX-11i) Nuova(E-8000)	

Descriptions	OEM Machines & Sensors	Price/USD
 OS-09 Response Time: <12s Warranty: 16months	Envitec Cells(OOM202-2) Maxtec(MAX-14) MSA Cells(806572) MSA Instruments(Mini-Ox IA) Nuova(E-09/ 0)	
 OS-10 Response Time: <12s Warranty: 16months	Analytical Industries(PSR-11-915-4) Envitec Cells(OOM110) GE Datex-Engström(Modulus SE7900) GE-Datex-Ohmeda Cells(6050-0004-110) GE-DatexOhmeda Instrum(Excel210-7900, Modulus SE1700, Smart Vent, S-5, Aespire/Avance, Airsys, Aestiva (3000)) Maxtec(MAX-10) Nuova(E-15/4) Teledyne Cells(R-30MED) Viamed(R-30V)	
 OS-11 Response Time: <12s Warranty: 16months	City Technologies(MOX-10) Envitec Cells(OOM202-2S) Nuova(E-110) Viamed(R-42V)	
 OS-12 Response Time: <12s Warranty: 16months	Analytical Industries(PSR-11-917) City Technologies(MOX-1) Envitec Cells(OOM102) GE-Datex-Ohmeda Cells(6051-0000-216) GE-DatexOhmeda Instrum(Centiva S / 5) Hamilton Cells(HM-12) Hamilton Instruments(Galileo, Raphael, Arabella, Multinex Analys) Hudson-RCI Cells(5566 , V-12) Maxtec(MAX-12A) Nuova(E-17/A2) Taema Instruments(Cesar (new)) Teledyne Cells(R-22MED) Viamed(R-22Med)	

Descriptions	OEM Machines & Sensors	Price/USD
 OS-13 Response Time: <15s Warranty: 12months	Analytical Industries(PSR-11-55) Hudson-RCI Cells(5500, V-05) Hudson-RCI Instruments (6477, 5577, 5584, 5590, 5595) Maxtec(MAX-9) MSA Cells(655265) Nuova(E-55/0) Teledyne Cells(R-13) Teledyne Instruments(#B 39636) Viamed(R-13)	
 OS-14 Response Time: <12s Warranty: 16months	Maxtec(MAX-13C) Nuova(E-14/0)	
 OS-15 Response Time: <12s Warranty: 16months	Analytical Industries(PSR-11-917J) City Technologies(MOX-2) Envitec Cells(OOM102-1) Maxtec(MAX-7) MSA Instruments(Mini-OX3000) Nuova(E-15/5) Taema Instruments(Alys, Clarys, YR010100) Teledyne Cells(R-17MED) Teledyne Instruments(TED AX 300) Viamed(R-26Med)	
 OS-16 Response Time: <12s Warranty: 16months	Maxtec(MAX-16) Nuova(E-16/0)	
 OS-16HT Response Time: <12s Warranty: 16months	Maxtec(MAX-16) Nuova(E-16/HTT)	

Descriptions		OEM Machines & Sensors	Price/USD
 OS-18T Response Time: <12s Warranty: 16months		MAXTEC(MAX-23)	
 OS-25 Response Time: <12s Warranty: 16months		Analytical Industries(PSR-11-917-J1) Maxtec(MAX-250ESF) Nuova(E-25/0)	
 OS-43 Response Time: <17s Warranty: 16months		GE-DatexOhmeda Instrum.(Giraffe) Maxtec(MAX-43) Nuova(E-43/0)	
 OS-47 Response Time: <12s Warranty: 12months		Maxtec(MAX-47) Nuova(E-47/0)	
 OS-48 Response Time: <12s Warranty: 12months		Analytical Industries(M-11-60) Maxtec(MAX-48) Nuova(E-60)	

Descriptions	OEM Machines & Sensors	Price/USD
 OS-80 Response Time: <12s Warranty: 12months	Envitec Cells(OOM204) Nuova(-17/A4)	

■ Oxygen Sensor Accessories

Product	Descriptions	Packaging	Price
 KLN-026 Warranty: 12months	15mm Tee-connector	1 pcs/bag	

1. The prices are ex-factory from Shenzhen, China.
2. Return for exchanging or refund is acceptable while items are found quality issue before use.