

ECG Cables and Leadwires

Features

- *Performance equivalent to OEMs*
- *Quality approved with FDA 510K and CE mark*
- *AHA and IEC codes available*
- *Multi-electrodes for best performance in the clinical environments*
- *Clear signals and secure connections to monitoring devices and patients*
- *6-12 months Warranty*



A comprehensive selection of replacement Cables and Leadwires available for these primary manufacturers



- | | |
|----------------------------|------------------------|
| ■ Dräger - Siemens | ■ Fufuda Danshi |
| ■ GE | ■ Medtronic |
| ■ Mindray Datascope | ■ Nihon Kohden |
| ■ Philips | ■ Spacelabs |

Comprehensive range of ECG Cables and Leadwires meeting the demands of all clinical purpose.

Shielded Cables and Leadwires

One-piece cable (2.5m+0.9m/1.2m)	Trunk Cables (2.5m)	Lead wires (0.9m/1.2m)
		

Din Safety Cables and Leadwires

Trunk cable (2.5m)	Lead with snap end (0.6m)	Leads w mini-grabber (0.6m)
		

Electrode Ends

Snap	Grabber	Mini Grabber	4mm Banana	3mm Needle
				

ECG Code/Lead Number

Lead#	AHA												IEC											
	RL	RA	LA	LL	V	V1	V2	V3	V4	V5	V6	N	R	L	F	C	C1	C2	C3	C4	C5	C6		
3		○	●	●									●	●	●									
4	●	○	●	●								●	●	●	●									
5	●	○	●	●	●							●	●	●	●	○								
6	●	○	●	●	●					●		●	●	●	●	○					●			
7	●	○	●	●	●				●	●		●	●	●	●	○				●	●			
10	●	○	●	●		●	●	●	●	●	●	●	●	●	●		●	●	●	●	●	●		

■ AAMI 6 Pin

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Trunk cable					
		3 lead used with Din Safety lead wire 2.5m	ET-3001-D-A	ET-3001-D-I	\$25.00
		5 lead used with Din Safety lead wire 2.5m	ET-5001-D-A	ET-5001-D-I	\$30.00
		3 lead used with Din Safety lead wire 2.5m	ET-3002-D-A	ET-3002-D-I	\$25.00
		5 lead used with Din Safety lead wire 2.5m	ET-5002-D-A	ET-5002-D-I	\$30.00
One-piece cable					
	Snap end	3 lead	EC-3001S-A	EC-3001S-I	\$25.00
	Grabber end	2.5m+0.9m	EC-3001G-A	EC-3001G-I	
	Snap end	5 lead	EC-5001S-A	EC-5001S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-5001G-A	EC-5001G-I	
	Snap end	3 lead	EC-3002S-A	EC-3002S-I	\$25.00
	Grabber end	2.5m+0.9m	EC-3002G-A	EC-3002G-I	
	Snap end	5 lead	EC-5002S-A	EC-5002S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-5002G-A	EC-5002G-I	






■ Artema

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
	Snap end	3 lead	EC-3010S-A	EC-3010S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3010G-A	EC-3010G-I	
	Snap end	5 lead	EC-5010S-A	EC-5010S-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5010G-A	EC-5010G-I	


■ Biolight

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
 	Snap end	3 lead 2.5m+0.9m	EC-3052RS-A	EC-3052RS-I	\$30.00
	Grabber end		EC-3052RG-A (15-027-0002)	EC-3052RG-I (15-027-0004)	
	Snap end	5 lead 2.5m+0.9m	EC-5052RS-A	EC-5052RS-I	\$35.00
	Grabber end		EC-5052RG-A (15-027-0001)	EC-5052RG-I (15-027-0003)	
 	Snap end	3 lead 2.5m+0.9m	EC-3053RS-A	EC-3053RS-I	\$30.00
	Grabber end		EC-3053RG-A (15-031-0014)	EC-3053RG-I (15-031-0013)	
	Snap end	3 lead 2.5m+0.9m	EC-5053RS-A	EC-5053RS-I	\$35.00
	Grabber end		EC-5053RG-A (15-031-0004)	EC-5053RG-I (15-031-0002)	

Bionet

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Lead wires					
	Snap end	3 lead wire	EL-BN3-90S-A	EL-BN3-90S-I	\$18.00
	Grabber end	0.9m	EL-BN3-90G-A	EL-BN3-90G-I	
	Snap end	5 lead wire	EL-BN5-90S-A	EL-BN5-90S-I	\$25.00
	Grabber end	0.9m	EL-BN5-90G-A	EL-BN5-90G-I	
Trunk cable					
		3 lead used with Bionet lead wire 2.5m	ET-3044R-BN-A	ET-3044R-BN-I	\$30.00
		5 lead used with Bionet lead wire 2.5m	ET-5044R-BN-A	ET-5044R-BN-I	\$35.00
		3 lead used with Bionet lead wire 2.5m	ET-3045-BN-A	ET-3045-BN-I	\$30.00
		5 lead used with Bionet lead wire 2.5m	ET-5045-BN-A	ET-5045-BN-I	\$35.00
One-piece cable					
	Snap end	3 lead	EC-3044RS-A	EC-3044RS-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3044RG-A	EC-3044RG-I	
	Snap end	5 lead	EC-5044RS-A	EC-5044RS-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5044RG-A	EC-5044RG-I	
	Snap end	3 lead	EC-3045S-A	EC-3045S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3045G-A	EC-3045G-I	
	Snap end	5 lead	EC-5045S-A	EC-5045S-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5045G-A	EC-5045G-I	

CasMed


Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
	Snap end	3 lead	EC-3047S-A	EC-3047S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3047G-A	EC-3047G-I	
	Snap end	3 lead	EC-3048S-A	EC-3048S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3048G-A	EC-3048G-I	

Orantech Inc.




T +86-755 2369 9939 F +86-755 2324 3594 E info@orantech.com W www.orantech.com

The difference you can sense

■ Chamcare

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
	Snap end	3 lead	EC-3064S-A	EC-3064S-I	\$25.00
	Grabber end	2.5m+0.9m	EC-3064G-A	EC-3064G-I	\$25.00

■ Colin-Omron

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)	
Lead wires-						
	Snap end	3 lead	EL-CL3-90S-A	EL-CL3-90S-I	\$18.00	
	Grabber end	0.9m	EL-CL3-90G-A	EL-CL3-90G-I		
	Snap end	5 lead	EL-CL5-90S-A	EL-CL5-90S-I	\$25.00	
	Grabber end	0.9m	EL-CL5-90G-A	EL-CL5-90G-I		
Trunk cable						
		3/5 lead used with Colin lead wire 2.5m	ET-5035-CL-A	ET-5035-CL-I	\$35.00	
		3/5 lead used with Colin lead wire 2.5m	ET-5036-CL-A	ET-5036-CL-I	\$35.00	
One-piece cable						
	Snap end	3 lead	EC-3035S-A	EC-3035S-I	\$30.00	
	Grabber end	2.5m+0.9m	EC-3035G-A	EC-3035G-I		
		Snap end	5 lead	EC-5035S-A	EC-5035S-I	\$35.00
		Grabber end	2.5m+0.9m	EC-5035G-A	EC-5035G-I	
		Snap end	3 lead	EC-3036S-A	EC-3036S-I	\$30.00
		Grabber end	2.5m+0.9m	EC-3036G-A	EC-3036G-I	
		Snap end	5 lead	EC-5036S-A	EC-5036S-I	\$35.00
		Grabber end	2.5m+0.9m	EC-5036G-A	EC-5036G-I	







Comen

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
	Snap end	3 lead	EC-3056S-A	EC-3056S-I	\$25.00
	Grabber end	2.5m+0.9m	EC-3056G-A	EC-3056G-I	
	Snap end	3 lead	EC-3057S-A	EC-3057S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3057G-A	EC-3057G-I	
	Snap end	3 lead	EC-5057S-A	EC-5057S-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5057G-A	EC-5057G-I	


Din safety




Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Lead wires-Din safety					
	Snap end	3 lead	EL-D3-60S-A	EL-D3-60S-I	\$6.00
	Grabber end	Din safety	EL-D3-60G-A	EL-D3-60G-I	
	Mini-Grabber	0.6m	EL-D3-60MG-A	EL-D3-60MG-I	\$6.00
	Snap end	4 lead	EL-D4-60S-A	EL-D4-60S-I	\$8.00
	Grabber end	Din safety 0.6m	EL-D4-60G-A	EL-D4-60G-I	\$8.00
	Snap end	5 lead	EL-D5-60S-A	EL-D5-60S-I	\$10.00
	Grabber end	Din safety	EL-D5-60G-A	EL-D5-60G-I	
	Mini-Grabber	0.6m	EL-D5-60MG-A	EL-D5-60MG-I	\$10.00
	Snap end	6 lead	EL-D6-60S-A	EL-D6-60S-I	\$12.00
	Grabber end	Din safety 0.6m	EL-D6-60G-A	EL-D6-60G-I	
	Snap end	7 lead	EL-D7-60S-A	EL-D7-60S-I	\$14.00
	Grabber end	Din safety 0.6m	EL-D7-60G-A	EL-D7-60G-I	
	Snap end	10 lead Din safety 0.6m	EL-D10-60S-A	EL-D10-60S-I	\$20.00
	Grabber end		EL-D10-60G-A	EL-D10-60G-I	

■ Draeger

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Lead wires-Dual Pin					
	Snap end	3 lead	EL-SM3-90S-A	EL-SM3-90S-I	\$18.00
	Grabber end	Dual-Pin 0.9m	EL-SM3-90G-A	EL-SM3-90G-I	
	Snap end	5 lead	EL-SM5-90S-A	EL-SM5-90S-I	\$25.00
	Grabber end	Dual-Pin 0.9m	EL-SM5-90G-A (5956466)	EL-SM5-90G-I	
	Snap end	3 lead	EL-SMB3-90S-A	EL-SMB3-90S-I	\$18.00
	Grabber end	Dual-Pin 0.9m	EL-SMB3-90G-A	EL-SMB3-90G-I	
Lead wires-Single Pin					
	Snap end	3 lead	EL-DG3-90S-A	EL-DG3-90S-I	\$25.00
	Grabber end	Single-Pin 0.9m	EL-DG3-90G-A (MS16231)	EL-DG3-90G-I (MS161159)	
	Snap end	5 lead	EL-DG5-90S-A	EL-DG5-90S-I	\$30.00
	Grabber end	Single-Pin 0.9m	EL-DG5-90G-A (MS16546)	EL-DG5-90G-I (MS16158)	
	Snap end	6 lead	EL-DG6-90S-A	EL-DG6-90S-I	\$35.00
	Grabber end	Single-Pin 0.9m	EL-DG6-90G-A (MS16547)	EL-DG6-90G-I (MS16157)	
One-piece cable					
	Snap end	3 lead	EC-3026S-A	EC-3026S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3026G-A	EC-3026G-I	
	Snap end	5 lead	EC-5026S-A	EC-5026S-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5026G-A	EC-5026G-I	


Draeger (continued)

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Trunk cable					
 	 	3 lead used with Din safety lead wire 2.5m	ET-3027-D-A	ET-3027-D-I	\$30.00
 	 	3 lead used with Dual Pin lead wire 2.5m	ET-3026-SM-A	ET-3026-SM-I	\$30.00
	 	5 lead used with Dual Pin lead wire 2.5m	ET-5026-SM-A	ET-5026-SM-I	\$35.00
 	 	3 lead used with Dual-Pin lead wire 2.5m	ET-3028-SM-A	ET-3028-SM-I	\$30.00
	 	5 lead used with Dual-Pin lead wire 2.5m	ET-5028-SM-A	ET-5028-SM-I	\$35.00

Multi-Parameter cable						
		MultiMed 5 Pod, for ECG, Resp., SpO2, and temp.,2.5m	ETM-SM-E5ST (3368391)	AHA/IEC are the same	\$200.00	
		MultiMed 5 Pod, for ECG, Resp., SpO2, and temp.,1.5m	ETM-SM15-E5ST (5950196)	AHA/IEC are the same	\$200.00	
		MultiMed puls, for ECG, Resp.,SpO2, and temp., 2.5m	ETM-DG-E6ST (MS20093)	AHA/IEC are the same	\$250.00	
		MultiMed puls OR, for ECG, Resp., SpO2, temp., 2.5m	ETM-DG-E6ST-OR (MS20094)	AHA/IEC are the same	\$300.00	

Extension Cable							
				6 lead wire 2.5m	ET-DG6-DG-A (MS16256)	ET-DG6-DG-I	\$40.00


■ ECAFIX

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
	Snap end	5 lead	EC-5065RS-A	EC-5065RS-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5065RG-A	EC-5065RG-I	







■ Emtel

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
	Snap end	3 lead	EC-3049RS-A	EC-3049RS-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3049RG-A	EC-3049RG-I	
	Snap end	5 lead	EC-5049RS-A	EC-5049RS-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5049RG-A	EC-5049RG-I	



■ Fukuda Denshi

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Trunk cable					
		3 lead used with Dual-Pin lead wire 2.5m	ET-3031R-SM-A (CI-700D-3)	ET-3031R-SM-I	\$30.00
		5 lead used with Dual-Pin lead wire 2.5m	ET-5031R-SM-A (CI-700D-5)	ET-5031R-SM-I	\$35.00
		3 lead used with Dual-Pin lead wire 2.5m	ET-3032-SM-A	ET-3032-SM-I	\$30.00
		5 lead used with Dual-Pin lead wire 2.5m	ET-5032-SM-A	ET-5032-SM-I	\$35.00
One-piece cable					
	Snap end	3 lead	EC-3031RS-A	EC-3031RS-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3031RG-A	EC-3031RG-I	
	Snap end	5 lead	EC-5031RS-A	EC-5031RS-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5031RG-A	EC-5031RG-I	
	Snap end	12 lead	EC-10031RS-A	EC-10031RS-I	\$60.00
	Grabber end	2.5m+0.9m	EC-10031RG-A	EC-10031RG-I	
	Snap end	3 lead	EC-3032S-A	EC-3032S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3032G-A	EC-3032G-I	
	Snap end	5 lead	EC-5032S-A	EC-5032S-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5032G-A	EC-5032G-I	
	Snap end	3 lead	EC-3076G-A	EC-3076G-I	\$50.00
	Grabber end	2.5m+0.9m	EC-3076S-A	EC-3076S-I	
	Snap end	5 lead	EC-5076S-A (CI0-07CTP-5NA)	EC-5076S-I	\$60.00
	Grabber end	2.5m+0.9m	EC-5076G-A	EC-5076G-I	






GE

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Lead wires-Multi-Link					
	Snap end	3 lead Multi-Link 0.9m	EL-MQ3-90S-A	EL-MQ3-90S-I (411203-003)	\$20.00
	Grabber end		EL-MQ3-90G-A (412682-001)	EL-MQ3-90G-I (412682-003)	
	Snap end	5 lead Multi-Link Limb 0.9m	EL-MQ5-90S-A (411202-003)	EL-MQ5-90S-I (411200-003)	\$25.00
	Grabber end		EL-MQ5-90G-A (412681-001)	EL-MQ5-90G-I (412681-003)	
	Snap end	5 lead Multi-Link Chest 0.9m	EL-MQV5-90S-A	EL-MQV5-90S-I	\$25.00
	Grabber end		EL-MQV5-90G-A (416467-001)	EL-MQV5-90G-I (414467-002)	
	Snap end	6 lead Multi-Link 0.9m	EL-MQS6-90S-A	EL-MQS6-90S-I (421931-001)	\$30.00
	Grabber end		EL-MQS6-90G-A	EL-MQS6-90G-I (421933-001)	
Lead wires-Datex					
	Snap end	3 lead 300 series 0.9m	EL-DX3-90S-A (545327)	EL-DX3-90S-I	\$20.00
	Grabber end		EL-DX3-90G-A (545317)	EL-DX3-90G-I (545315)	
	Snap end	5 lead 300 series 0.9m	EL-DX5-90S-A	EL-DX5-90S-I	\$25.00
	Grabber end		EL-DX5-90G-A	EL-DX5-90G-I (545316)	

GE (continued)

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Trunk cable-Marquette					
		3/5 lead used with Multi-Link lead wire 2.5m	ET-5004-MQ-A (2017003-001)	ET-5004-MQ-I (2017003-003)	\$35.00
		3 lead used with Din Safety lead wire 2.5m	ET-3004-D-A (2017004-001)	ET-3004-D-I (2017004-003)	\$30.00
		6 lead used with Multi-Link lead wire 2.5m	ET-6004-MQ-A (2017005-001)	ET-6004-MQ-I (2017005-003)	\$40.00
		12 lead used with Multi-Link lead wire 2.5m	ET-10004-MQ-A (2017006-002)	ET-10004-MQ-I (2017006-003)	\$50.00
One-piece cable-Marquette					
	Snap end	3 lead 2.5m+0.9m	EC-3004S-A (2106309-001)	EC-3004S-I	\$30.00
	Grabber end		EC-3004G-A (2021141-001)	EC-3004G-I (2021141-002)	\$30.00
	Snap end	5 lead 2.5m+0.9m	EC-5004S-A	EC-5004S-I	\$35.00
	Grabber end		EC-5004G-A	EC-5004G-I	
	Snap end	10lead 2.5+0.9/1.2m	EC-10004S-A	EC-10004S-I	\$60.00
	Grabber end		EC-10004G-A	EC-10004G-I	


GE (continued)

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Trunk cable-Datex					
		3 lead used with 300 series lead wire 2.5m	ET-3008R-DX-A (545302)	ET-3008R-DX-I (545300)	\$30.00
		5 lead used with 300 series lead wire 2.5m	ET-5008R-DX-A (545303)	ET-5008R-DX-I (545301)	\$35.00
		3 lead used with 300 series lead wire 2.5m	ET-3009-DX-A	ET-3009-DX-I	\$30.00
		5 lead used with 300 series lead wire 2.5m	ET-5009-DX-A	ET-5009-DX-I	\$35.00
One-piece cable-Datex					
	Snap end	3 lead	EC-3008RS-A	EC-3008RS-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3008RG-A	EC-3008RG-I	
	Snap end	5 lead	EC-5008RS-A	EC-5008RS-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5008RG-A	EC-5008RG-I	
	Snap end	3 lead	EC-3009S-A	EC-3009S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3009G-A	EC-3009G-I	
	Snap end	5 lead	EC-5009S-A	EC-5009S-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5009G-A	EC-5009G-I	
One-piece cable-Hellige					
	Snap end	3 lead	EC-3021RS-A	EC-3021RS-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3021RG-A	EC-3021RG-I	
	Snap end	5 lead	EC-5021RS-A	EC-5021RS-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5021RG-A	EC-5021RG-I	

Huntleigh

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
	Snap end	3 lead	EC-3067S-A	EC-3067S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3067G-A	EC-3067G-I	



■ Infinum

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
	Snap end	3 lead	EC-3058RS-A	EC-3058RS-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3058RG-A	EC-3058RG-I	
	Snap end	5 lead	EC-5058RS-A	EC-5058RS-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5058RG-A	EC-5058RG-I	






■ Kontron

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
	Snap end	3 lead	EC-3029RS-A	EC-3029RS-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3029RG-A	EC-3029RG-I	
	Snap end	5 lead	EC-5029RS-A	EC-5029RS-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5029RG-A	EC-5029RG-I	
	Snap end	3 lead	EC-3030S-A	EC-3030S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3030G-A	EC-3030G-I	


■ Medtronic-Physio Control

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Lead wires					
	Snap end	6 lead wire Medtronic	EL-MDV6-90S-A (11110-000022)	EL-MDV6-90S-I	\$35.00
	Grabber end	0.9m	EL-MDV6-90G-A	EL-MDV6-90G-I	\$35.00
One-piece cable					
	Snap end	3 lead 2.5m+0.9m	EC-3039RS-A (11110-000030)	EC-3039RS-I	\$30.00
	Snap end	4 lead 2.5m+0.9m	EC-4039RS-A (11110-000018)	EC-4039RS-I	\$60.00
	Snap end	5 lead 2.5m+0.9m	EC-5039RS-A (11110-000067)	EC-5039RS-I	\$35.00
	Snap end	12 lead 2.5m+0.9m	EC-10039RS-A	EC-10039RS-I	\$60.00





Mindray>Datascope

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Lead wires-Mindray					
	Snap end	3 lead used with	EL-MR3-90S-A	EL-MR3-90S-I	\$18.00
	Grabber end	MR lead wire,0.9m	EL-MR3-90G-A	EL-MR3-90G-I	
	Mini-grabber end	3 lead used with MR lead wire,0.9m	EL-MR3-90MG-A	EL-MR3-90MG-I	\$20.00
	Snap end	5 lead used with MR lead wire	EL-MR5-90S-A (EL6501B)	EL-MR5-90S-I	\$25.00
	Grabber end	0.9m	EL-MR5-90G-A (EL6501A)	EL-MR5-90G-I	
	Snap end	3 lead used with MRB lead wire	EL-MRB3-90S-A (EL6301B)	EL-MRB3-90S-I	\$20.00
	Grabber end	0.9m	EL-MRB3-90G-A (EL6301A)	EL-MRB3-90G-I	
	Snap end	5 lead used with MRA lead wire	EL-MRA5-90S-A	EL-MRA5-90S-I	\$30.00
	Grabber end	Limb 0.9m	EL-MRA5-90G-A (0010-30-42906)	EL-MRA5-90G-I (0010-30-42902)	
	Snap end	5 lead used with	EL-MRV5-90S-A	EL-MRV5-90S-I	\$30.00
	Grabber end	MRV lead wire Chest, 0.9m	EL-MRV5-90P-A (0010-30-42908)	EL-MRV5-90P-I (0010-30-42904)	
Lead wires-Datascope					
	Snap end	3 lead used with DT lead wire	EL-DT3-90S-A (0012-00-1503-05)	EL-DT3-90S-I	\$20.00
	Grabber end	0.9m	EL-DT3-90G-A	EL-DT3-90G-I	
	Snap end	5 lead used with DT lead wire	EL-DT5-90S-A (0012-00-1503-00)	EL-DT5-90S-I	\$25.00
	Grabber end	0.9m	EL-DT5-90G-A	EL-DT5-90G-I	
Trunk cable-Datascope					
		5 lead used with DT lead wire,2.5m	ET-5018-DT-A	ET-5018-DT-I	\$35.00
One-piece cable-Datascope					
	Snap end	3 lead	EC-3017RS-A	EC-3017RS-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3017RG-A	EC-3017RG-I	
	Snap end	3 lead	EC-3018RS-A	EC-3018RS-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3018RG-A	EC-3018RG-I	
	Snap end	5 lead	EC-5018RS-A	EC-5018RS-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5018RG-A	EC-5018RG-I	


■ Mindray>Datascope (continued)

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Trunk cable-Mindray					
 		3 lead used with MR lead wire 2.5m	ET-3001-MR-A (EV6130N)	ET-3001-MR-I	\$25.00
			ET-3001R-MR-A (EV6102)	ET-3001R-MR-I	
 		3 lead used with MR lead wire 2.5m	ET-3016-MR-A (EV6230N)	ET-3016-MR-I	\$30.00
			ET-3016R-MR-A (EV6202)	ET-3016R-MR-I	
		5 lead used with MR lead wire 2.5m	ET-5016-MR-A (EV6250N)	ET-5016-MR-I	\$35.00
			ET-5016R-MR-A (EV6201)	ET-5016R-MR-I	\$35.00
		12 lead used with MRA lead wire, 2.5m	ET-10016R-MRA-A (0010-30-42721)	ET-10016R-MRA-I	\$50.00
One-piece cable-Mindray					
 	Snap end	3 lead 2.5m+0.9m	EC-3016S-A (EA231B)	EC-3016S-I	\$30.00
	Grabber end		EC-3016G-A (EA231A)	EC-3016G-I	
	Snap end	5 lead 2.5m+0.9m	EC-5016S-A (EA6251B)	EC-5016S-I	\$35.00
	Grabber end		EC-5016G-A (EA6251A)	EC-5016G-I	
	Snap end	10lead 2.5+0.9/1.2m	EC-10016G-A	EC-10016G-I	\$50.00
	Grabber end		EC-10016S-A	EC-10016S-I	
	Snap end	12 lead 2.5m+0.9m	EC-10016RS-A	EC-10016RS-I	\$50.00
	Grabber end		EC-10016RG-A	EC-10016RG-I	








■ Mennen

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Lead wires					
	Snap end	3 lead wire	EL-MN3-90S-A	EL-MN3-90S-I	\$20.00
	Grabber end	Mennen 0.9m	EL-MN3-90G-A	EL-MN3-90G-I	
	Snap end	5 lead wire	EL-MN5-90S-A	EL-MN5-90S-I	\$25.00
	Grabber end	Mennen 0.9m	EL-MN5-90G-A	EL-MN5-90G-I	
Trunk Cable					
		3 lead used with Mennen lead wire 2.5m	ET-3051R-MN-A	ET-3051R-MN-I	\$25.00
		5 lead used with Mennen lead wire 2.5m	ET-5051R-MN-A	ET-5051R-MN-I	\$30.00
One-piece Cable					
	Snap end	3 lead	EC-3051RS-A	EC-3051RS-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3051RG-A	EC-3051RG-I	
	Snap end	5 lead	EC-5051RS-A	EC-5051RS-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5051RG-A	EC-5051RG-I	




■ MEK

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
	Snap end	3 lead	EC-3033S-A	EC-3033S-I	\$25.00
	Grabber end	2.5m+0.9m	EC-3033G-A	EC-3033G-I	
	Snap end	5 lead	EC-5033S-A	EC-5033S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-5033G-A	EC-5033G-I	
	Snap end	3 lead	EC-3034S-A	EC-3034S-I	\$25.00
	Grabber end	2.5m+0.9m	EC-3034G-A	EC-3034G-I	
	Snap end	5 lead	EC-5034S-A	EC-5034S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-5034G-A	EC-5034G-I	






Nihon Kohden

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Lead wires-Nihon Kohden					
	Snap end	3 lead NKA 0.9m	EL-NKA3-90S-A (BR-019PA)	EL-NKA3-90S-I (BR-019P)	\$20.00
	Grabber end		EL-NKA3-90G-A (BR-018PA)	EL-NKA3-90G-I (BR-018P)	
	Snap end	5 lead NKA 0.9m	EL-NKA5-90S-A (BR-021PA)	EL-NKA5-90S-I (BR-021P)	\$25.00
	Grabber end		EL-NKA5-90G-A (BR-020PA)	EL-NKA5-90G-I (BR-020P)	
	Snap end	3 lead NKB 0.9m	EL-NKB3-90S-A (BR-913PA)	EL-NKB3-90S-I (BR-913P)	\$20.00
	Grabber end		EL-NKB3-90G-A (BR-903PA)	EL-NKB3-90G-I (BR-903P)	
	Snap end	6 lead NKB 0.9m	EL-NKB6-90S-A (BR-916PA)	EL-NKB6-90S-I (BR-916P)	\$25.00
	Grabber end		EL-NKB6-90G-A (BR-906PA)	EL-NKB6-90G-I (BR-906P)	
Trunk cable					
		3/5 lead used with NKA lead wire 2.5m	ET-5011R-NKA-A	ET-5011R-NKA-I	\$35.00
		3/5 lead used with NKA lead wire 2.5m	ET-5012R-NKA-A	ET-5012R-NKA-I	\$35.00
		3/5 lead used with NKB lead wire 2.5m	ET-3015R-NKB-A (JC-103TA)	ET-3015R-NKB-I (JC-103T)	\$30.00
		3/6 lead used with NKB lead wire 2.5m	ET-6013-NKB-A (JC-906PA)	ET-6013-NKB-I (JC-906P)	\$35.00
		3 lead used with NKB lead wire 2.5m	ET-3014R-NKB-A	ET-3014R-NKB-I	\$35.00



■ Nihon Kohden (continued)

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)		
One-piece cable							
		Snap end	3 lead	EC-3011RS-A	EC-3011RS-I	\$30.00	
		Grabber end	2.5m+0.9m	EC-3011RG-A	EC-3011RG-I		
			Snap end	5 lead	EC-5011RS-A	EC-5011RS-I	\$35.00
			Grabber end	2.5m+0.9m	EC-5011RG-A	EC-5011RG-I	
			Snap end	3 lead	EC-3012RS-A	EC-3012RS-I	\$30.00
			Grabber end	2.5m+0.9m	EC-3012RG-A	EC-3012RG-I	
			Snap end	5 lead	EC-5012RS-A	EC-5012RS-I	\$35.00
			Grabber end	2.5m+0.9m	EC-5012RG-A	EC-5012RG-I	
			Snap end	3 lead	EC-3014RS-A (BC-763VA)	EC-3014RS-I (BC-763V)	\$30.00
			Grabber end	2.5m+0.9m	EC-3014RG-A (BC-753VA)	EC-3014RG-I (BC-753V)	\$30.00
		Snap end	12 lead	EC-10013S-A	EC-10013S-I	\$60.00	
		Grabber end	2.5m+0.9m	EC-10013G-A (BJ-900PA)	EC-10013G-I (BJ-900P)		
		Snap end	5lead	EC-5014RG-A	EC-5014RG-I	\$35.00	
		Grabber end	2.5m+0.9m	EC-5014RS-A	EC-5014RS-I	\$35.00	

Philips

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Lead wires-Twin Pin					
	Snap end	3 lead Twin-Pin 0.9m	EL-AA3-90S-A (M1605A)	EL-AA3-90S-I (M1615A)	\$20.00
	Grabber end		EL-AA3-90G-A (M1603A)	EL-AA3-90G-I (M1613A)	
	Snap end	5 lead Twin-Pin 0.9m	EL-AA5-90S-A (M1625A)	EL-AA5-90S-I (M1635A)	\$25.00
	Grabber end		EL-AA5-90G-A (M1623A)	EL-AA5-90G-I (M1633A)	
Lead wires-IntelliVue					
	Snap end	3 lead IntelliVue 0.9m	EL-PH3-90S-A (M1673A)	EL-PH3-90S-I (M1674A)	\$20.00
	Grabber end		EL-PH3-90G-A (M1671A)	EL-PH3-90G-I (M1672A)	
	Snap end	5 lead IntelliVue Limb 0.9m	EL-PH5-90S-A (M1644A)	EL-PH5-90S-I (M1645A)	\$25.00
	Grabber end		EL-PH5-90G-A (M1968A)	EL-PH5-90G-I (M1971A)	
	Snap end	5 lead IntelliVue Chest 0.9m	EL-PHV5-90S-A (M1602A)	EL-PHV5-90S-I (M1604A)	\$25.00
	Grabber end		EL-PHV5-90G-A (M1976A)	EL-PHV5-90G-I (M1978A)	

Philips (continued)


Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Trunk cable					
		3 lead used with Twin-Pin lead wire 2.5m	ET-3006R-AA-A	ET-3006R-AA-I (M1590A)	\$30.00
		5 lead used with Twin-Pin lead wire 2.5m	ET-5006R-AA-A	ET-5006R-AA-I	\$35.00
One-piece cable					
	Snap end	3 lead 2.5m+0.9m	EC-3006RS-A (M1733A)	EC-3006RS-I (M1735A)	\$30.00
	Snap end	5 lead 2.5m+0.9m	EC-5006RS-A (M1734A)	EC-5006RS-I (M1736A)	\$35.00
Trunk cable					
		3 lead used with Twin-Pin lead wire 2.5m	ET-3005R-AA-A (M1500A)	ET-3005R-AA-I (M1510A)	\$30.00
		3 lead used with IntelliVue lead wire 2.5m	ET-3005R-PH-A (M1669A)	ET-3005R-PH-I (M1669A)	\$30.00
		5 lead used with Twin-Pin lead wire 2.5m	ET-5005R-AA-A (M1520A)	ET-5005R-AA-I (M1530A)	\$35.00
		5 lead used with IntelliVue lead wire 2.5m	ET-5005R-PH-A (M1668A)	ET-5005R-PH-I (M1668A)	\$35.00
		12 lead used with IntelliVue lead wire 2.5m	ET-10005R-PH-A (M1949A, M1663A)	ET-10005R-PH-I (M1949A, M1663A)	\$50.00
One-piece cable					
	Snap end	3 lead 2.5m+0.9m	EC-3005RS-A	EC-3005RS-I	\$30.00
	Grabber end		EC-3005RG-A (989803143181)	EC-3005RG-I (989803143171)	
	Snap end	5 lead 2.5m+0.9m	EC-5005RS-A	EC-5005RS-I	\$35.00
	Grabber end		EC-5005RG-A	EC-5005RG-I (989803143191)	
	Snap end	10 lead 2.5+0.9/1.2m	EC-10005RS-A	EC-10005RS-I	\$60.00
	Grabber end		EC-10005RG-A	EC-10005RG-I	

Orantech Inc.








T +86-755 2369 9939 F +86-755 2324 3594 E info@orantech.com W www.orantech.com

The difference you can sense



■ Schiller

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
	Snap end	3 lead	EC-3040S-A	EC-3040S-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3040G-A	EC-3040G-I	
	Snap end	3 lead	EC-3041RS-A	EC-3041RS-I	\$30.00
	Grabber end	2.5m+0.9m	EC-3041RG-A	EC-3041RG-I	
	Snap end	5 lead	EC-5075RS-A	EC-5075RS-I	\$35.00
	Grabber end	2.5m+0.9m	EC-5075RG-A	EC-5075RG-I	

■ Spacelabs

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Lead wires-TruLink					
	Snap end	3 lead	EL-SL3-90S-A (700-0007-01/03)	EL-SL3-90S-I	\$18.00
	Grabber end	TruLink 0.9m	EL-SL3-90G-A	EL-SL3-90G-I (700-0006-01/03)	
	Snap end	5 lead	EL-SL5-90S-A 700-0006-09/11)	EL-SL5-90S-I	\$25.00
	Grabber end	TruLink 0.9m	EL-SL5-90G-A	EL-SL5-90G-I	
Trunk Cable					
		3 lead used with TruLink lead wire 2.5m	ET-3007R-SL-A (700-0008-08)	ET-3007R-SL-I (700-0008-09)	\$30.00
		5 lead used with TruLink lead wire 2.5m	ET-5007R-SL-A (700-0008-06)	ET-5007R-SL-I (700-0008-07)	\$35.00
		3 lead used with Safety Din lead wire 0.6m	ET-3007R-D-A (700-0008-10)	ET-3007R-D-I	\$30.00
One-piece Cable					
	Snap end	12 lead	EC-10007RS-A	EC-10007RS-I	\$60.00
	Grabber end	2.5m+0.9m	EC-10007RG-A	EC-10007RG-I	


Siemens CT

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
Trunk Cables					
		10 lead 2.5m	ET-10046-SMCT-A (6644541)	ET-10046-SMCT-I	\$50.00

Welch Allyn

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
	Snap end	3 lead 2.5m+0.9m	EC-3019RS-A (008-0315-00)	EC-3019RS-I	\$30.00
	Grabber end		EC-3019RG-A	EC-3019RG-I	
	Snap end	3 lead 2.5m+0.9m	EC-3020RS-A (008-0880-00)	EC-3020RS-I (008-0880-01)	\$30.00
	Grabber end		EC-3020RG-A	EC-3020RG-I	
	Snap end	5 lead 2.5m+0.9m	EC-5020RS-A (008-0879-00)	EC-5020RS-I	\$35.00
	Grabber end		EC-5020RG-A	EC-5020RG-I	

Zoll

Proximal end Connector	Distal end Connector	Lead number	Part Number/AHA (OEM PN)	Part Number/IEC (OEM PN)	Price (USD)
One-piece cable					
	Snap end	12 lead 2.5m+0.9m	EC-10037RS-A	EC-10037RS-I	\$60.00
	Grabber end		EC-10037RG-A	EC-10037RG-I	
	Snap end	12 lead 2.5m+0.9m	EC-10038RS-A (103801)	EC-10038RS-I (103802)	\$60.00
	Grabber end		EC-10038RG-A	EC-10038RG-I	